

### 3.3 Sample Collection, Storage and Preservation

Figure 3 details required containers, sample volumes, preservation techniques, and holding times for proper sample collection. A discussion is included pertaining to sample collection, preservation and storage before analysis.

#### Sample Collection

RMBEL supplies clients with various field services and sample collections. These field services include, but are not limited to, Secchi disk readings, lake monitoring sampling, ground water monitoring, surface water monitoring and private drinking water collections. All personnel are trained with sampling protocol and approved standard operating procedures to ensure correct sample handling. If a situation arises where the sampling deviates from said protocols or SOPs, the deviation is recorded in detail on the sampling form and chain-of-custody form. The client is notified of the deviation if applicable.

**Appendix D** lists the Field Services Standard Operating Procedures. **Appendix E** lists the Biological Monitoring Standard Operating Procedures.

#### Holding Times

The integrity of the sample is dependent on strict adherence to standard operating procedures in place in the laboratory. RMBEL maintains and recognizes the importance of sample holding times. Sample holding times are determined by the analytical method. Holding times are noted in Figure 3. Samples received after the holding time has expired are analyzed only under direct request from the client. Appropriate sample and/or data qualifiers are noted on the final report, as applicable (**Appendix J**).

#### Preservation Techniques

Proper preservation practices must be followed and can be found for each analyte in Figure 3. Following collection and during transportation, samples should be kept at 6°C or on ice. Samples requiring preservation should be preserved as soon as possible after collection to maintain the integrity of the sample. Complete collection, storage, and preservation guidelines can be found in 40CFR Part 136.

NOTE: Complete and certain preservation of samples, regardless of source, is a practical impossibility. Regardless of the sample nature, complete stability for every constituent can never be fully attained. At best, sample preservation only slows the biological and chemical changes that inevitably continue after the sample is collected.

Methods of preservation are intended to retard biological action, retard hydrolysis of chemical compounds and complexes, and reduce volatility of constituents. Preservation methods are limited to pH control, chemical addition, amber or opaque bottles, filtration, refrigeration, and freezing.

To minimize the potential for volatilization or biodegradation between sampling and analysis, keep the sample as cool as possible without freezing. Preferably, pack samples in crushed or cubed ice, or a commercial ice substitute before shipment. Avoid using dry ice because it may freeze the samples and cause the containers to break. Dry ice may also affect the pH of the samples. Keep composite samples cool with ice or a refrigeration system set at 6°C during collection. Analyze the samples as quickly as possible upon arrival at the laboratory. If immediate analysis is not possible, storage at 2°-4° C is recommended for most samples.

Use chemical preservation only when it is shown not to interfere with the method of analysis. When chemicals are used, add them to the sample bottle so that all sample portions are evenly preserved as soon as collected. No single preservation method is entirely satisfactory; choose the preservation with regard to the analyses being made. Because a preservation method for one analysis may interfere with the preservation for another, samples for multiple determinations may need to be split and preserved separately. All methods of preservation may be inadequate when applied to suspended matter. DO NOT use formaldehyde as a preservative because it affects numerous analyses.

The previous discussion is by no means exhaustive and/or comprehensive. Clearly, it is impossible to prescribe absolute rules for preventing all possible changes. Additional advice can be found in the analytical methods and/or SOPs under

# COPY

individual determination, but to a large degree the dependability of an analytical determination rests on the experience and good judgment of the person collecting the sample.

### 3.4 Sample Storage and Disposal

#### Sample Storage Prior to Analysis

All samples are properly stored from the time they arrive at the laboratory to disposal. Samples are refrigerated at 4°C prior to analysis unless method SOPs indicate other storage conditions.

**Nature of Sample Changes:** Some analyses are more likely than others to be affected by storage before analysis. Certain cations are subject to loss by absorption or ion exchange with the walls of glass containers. These include aluminum, cadmium, chromium, copper, iron, lead, manganese, silver, and zinc, which are best collected in a separate clean bottle and acidified with nitric acid to a pH below 2.0. This minimizes precipitation and absorption to container walls.

Temperature, pH, and dissolved oxygen are best determined in the field. Temperature changes quickly and pH may change significantly in a matter of minutes. Dissolved gases (oxygen, carbon dioxide) may be lost very quickly. With the changes in the pH-alkalinity-carbon dioxide balance, calcium carbonate may precipitate and cause a decrease in the values for calcium and total hardness.

Iron and manganese are readily soluble in their lower oxidation states but are relatively insoluble in their higher oxidation states, therefore, these cations may precipitate or they may dissolve from sediment depending on the redox potential of the sample. Microbial action may be responsible for the changes in the nitrate-nitrite-ammonia content, or decreases in phenol concentration. It may also be responsible for reducing sulfate to sulfide in Biological Oxygen Demand (BOD). Residual chlorine is reduced to chloride. Sulfide, sulfite, ferrous iron, iodide, and cyanide may be lost through oxidation. Color, Odor, and Turbidity may increase, decrease, or change quality. Sodium, silica, and boron may be leached from a glass container. Hexavalent chromium may be reduced to chromic acid.

Biological changes taking place in a sample may change the oxidation state of some constituents. Soluble constituents may be converted to organic bound material in cell structures, or cell lysis may result in release of cellular material into solution. The well-known nitrogen and phosphorus cycles are examples of biological influence on sample composition.

Zero headspace is important in preservation of samples with volatile organics. Avoid loss of volatile materials by filling sample containers completely; achieve this by filling the bottle until there is a positive meniscus before capping and sealing, being careful not to rinse out the acid preservation or come into contact with the acidified sample. Serum vials with septum caps are particularly useful in that a sample portion for analysis can be taken through the cap with a syringe.

In general, the shorter the time that elapses between the collection of a sample and its analysis, the more reliable the analytical results. For certain constituents and physical values, immediate analysis in the field is required. For samples composited in the field, use the time at the end of the composition collection as the sample collection time.

It is impossible to state exactly how much time may be allowed between sample collection and analysis. Changes occurring in the sample depend on the character of the sample, the analysis to be made, and the conditions of storage. Changes caused by the growth of microorganisms are greatly retarded by keeping the sample in the dark and at a low temperature. When the interval between sample collection and analysis is long enough to produce changes in either the concentration or the physical state of the constituent to be measured, follow the preservation practices stated in Figure 3.

#### Sample Storage

Following sample receipt and proper preparations by the laboratory, all samples are stored separate from standards or reagents. A temperature of 4° C is maintained and monitored 24 hours a day by digital controls. AM and PM temperature readings are maintained in a daily temperature log. Samples are kept in the cooler until time of analysis. Samples that require analysis at room temperature (25° C) are brought out of refrigeration and allowed to warm to the desired temperature before beginning analysis. When all analyses are performed on a sample, it is transferred to a cooler for cold storage until the time of disposal.

# COPY

## Sample Disposal

All samples are disposed in an environmentally sound manner or returned to the client upon request. Samples may also be returned to the client if they contain hazardous wastes. Samples that contain solvents are evaporated in a fume hood unless the vapors create an environmental concern. Samples that contain hazardous materials, such as mercury and silver, are sent to a hazardous waste collection facility that has a Very Small Quantity Generator (VSQG) program. The remaining samples are stored in a holding tank that is transferred to a Wastewater Treatment Facility for proper disposal. The pH of the waste is verified and adjusted to be above pH 2.0 before being added to the holding tank. Refer to RMBEL SOP 03-005 for Waste Disposal procedures. Samples are typically stored for 30 days following the report date unless otherwise instructed by the client.

**Figure 3: Sampling Requirements**

Parameter	Matrix	Container	Minimum Sample Size	Preservative	Holding Time
<b>Microbiology</b>					
Coliform, Total or Fecal	Water & Wastewater	Sterile Container	100 mL	Refrigerate 6 C	8 Hrs EPA compliance 24 Hours
Coliform, Total or Fecal, Chlorinated Water	Water & Wastewater	Sterile Container w/ thiosulfate	100 mL	Refrigerate 6 C	8 Hrs EPA compliance 24 Hours
E. coli	Water & Wastewater	Sterile Container	100 mL	Refrigerate 6 C	8 Hrs EPA compliance 24 Hours
Total Coliform or E. coli	Drinking Water	Sterile Container	100 mL	Refrigerate 6 C	8 Hrs EPA compliance 30 hours
<b>Inorganics and Nutrients</b>					
Acidity, Free	Water & Wastewater	Plastic or Glass	500 mL	Cool 6 C	14 Days
Alkalinity, Total	Water & Wastewater	Plastic or Glass	500 mL	Cool 6 C	14 Days
BOD	Water & Wastewater	Plastic or Glass	650 mL	Cool 6 C	48 Hours
BOD, Carbonaceous	Water & Wastewater	Plastic or Glass	650 mL	Cool 6 C	48 Hours
Chloride	Water & Wastewater	Plastic or Glass	100 mL	None	28 Days
	Solid	Plastic	50 mL	None	28 Days
Chlorine, Residual	Water & Wastewater	Plastic or Glass	500 mL	None	Analyze Immediately
Chlorophyll-a	Water	Amber Glass	1000 mL	Cool 6 C, Dark	48 Hours
Chromium, Hexavalent	Water & Wastewater	Plastic or Glass, rinsed w/ 1:1 HNO <sub>3</sub>	200 mL	Cool 4 C	24 Hours
COD	Water & Wastewater	Plastic or Glass	100 ml	H <sub>2</sub> SO <sub>4</sub> to pH < 2, Cool 6C	28 Days
	Solid	Plastic	50 mL	None	28 Days
Color	Water & Wastewater	Plastic or Glass	250 mL	Cool 6 C	48 Hours
Conductance, Specific	Water & Wastewater	Plastic or Glass	500 ml	Cool 6 C	28 Days
	Solid	Plastic	50 mL	None	28 Days
Cyanide, Total	Water & Wastewater	Plastic or Glass	500 mL	NaOH to pH > 12, Cool 6C	14 Days
	Solid	Plastic	50 mL	None	14 Days
Fluoride	Water & Wastewater	Plastic	100 mL	None, Cool 6C	28 Days
	Solid	Plastic	50 mL	None	28 Days
Hardness, Total	Water & Wastewater	Plastic or Glass	500 mL	HNO <sub>3</sub> to pH < 2	6 Months
Iodide	Water & Wastewater	Plastic or Glass	500 mL	Cool 6 C	24 Hours
Nitrogen, Kjeldahl Total	Water & Wastewater	Plastic or Glass	500 mL	H <sub>2</sub> SO <sub>4</sub> to pH < 2, Cool 6C	28 Days
	Solid	Plastic	50 mL	Cool 6 C	28 Days
Nitrogen, Nitrate + Nitrite	Water & Wastewater	Plastic or Glass	250 mL	H <sub>2</sub> SO <sub>4</sub> to pH < 2, Cool 6C	28 Days
	Solid	Plastic	50 mL	Cool 6 C	28 Days
Nitrogen, Nitrate	Water & Wastewater	Plastic or Glass	250 mL	Cool 6 C	48 Hours
	Solid	Plastic	50 mL	Cool 6 C	48 Hours

# COPY

Figure 3: continued

Parameter	Matrix	Container	Minimum Sample Size	Preservative	Holding Time
Nitrogen, Nitrite	Water & Wastewater	Plastic or Glass	250 mL	Cool 6 C	48 Hours
	Solid	Plastic	50 mL	Cool 6 C	48 Hours
Nitrogen, Total	Water and Wastewater	Plastic or Glass	250 mL	H <sub>2</sub> SO <sub>4</sub> to pH < 2, Cool 6C	
Nitrogen, Ammonia	Water & Wastewater	Plastic or Glass	500 mL	H <sub>2</sub> SO <sub>4</sub> to pH < 2, Cool 6C	28 Days
	Solid	Plastic	50 mL	Cool 6 C	28 Days
Oil & Grease	Water & Wastewater	Amber Glass	1000 mL	H <sub>2</sub> SO <sub>4</sub> to pH < 2 Cool 6 C	28 Days
Organic Carbon, Total	Water & Wastewater	Plastic or Glass	500 mL	H <sub>2</sub> SO <sub>4</sub> to pH < 2, Cool 6C	28 Days
	Solid	Plastic	50 mL	Cool 6 C	28 Days
Oxygen, Dissolved	Water & Wastewater	Glass bottle w/ top, no headspace	500 mL	None	Analyze Immediately
pH	Water & Wastewater	Plastic or Glass	100 mL	None	Analyze in 15min
	Solid	Plastic	50 mL	None	Analyze in 15min
Phaeophytin	Water	Amber Glass	1000 mL	Cool 6 C	48 Hours
Phenolics, Total	Water & Wastewater	Amber, Plastic, or Glass	1000 mL	H <sub>2</sub> SO <sub>4</sub> to pH <2, Cool 6C	28 Days
Phosphorus, Total	Water and Wastewater	Plastic or Glass	500 mL	H <sub>2</sub> SO <sub>4</sub> to pH <2, Cool 6C	28 Days
	Solid	Plastic	50 mL	Cool 6 C	28 Days
Phosphorus, Ortho	Water & Wastewater	Plastic or Glass	500 mL	Filter, Cool 6 C	48 Hours, Filter within 15 minutes
	Solid	Plastic	50 mL	None	48 Hours
Salinity	Water & Wastewater	Glass	500 mL	None	Analyze Immediately
Settleable Matter	Water & Wastewater	Plastic or Glass	1000 mL	Cool 6 C	48 Hours
Silica	Water & Wastewater	Plastic	500 mL	Cool 6 C	28 Days
Solids, Total	Water & Wastewater	Plastic or Glass	500 mL	Cool 6 C	7 Days
	Solid	Plastic	50 mL	None	7 Days
Solids, Total Volatile	Water & Wastewater	Plastic or Glass	500 mL	Cool 6 C	7 Days
	Solid	Plastic	50 mL	None	7 Days
Solids, Total Suspended	Water & Wastewater	Plastic or Glass	500 mL	Cool 6 C	7 Days
Solids, Total Settleable	Water & Wastewater	Plastic or Glass	500 mL	Cool 6 C	7 Days
Solids, Total Dissolved	Water & Wastewater	Plastic or Glass	500 mL	Cool 6 C	7 Days
Solids, Suspended Volatile	Water & Wastewater	Plastic or Glass	500 mL	Cool 6 C	7 Days
Sulfate	Water & Wastewater	Plastic or Glass	500 mL	Cool 6 C	28 Days
	Solid	Plastic	50 mL	None	28 Days
Sulfide	Water & Wastewater	Amber Glass	500 mL	Cool 6 C 2 mL Zinc Acetate	7 Days
	Solid	Plastic	50 mL	None	7 Days
Sulfite	Water & Wastewater	Plastic or Glass	500 mL	None	Analyze Immediately
Surfactants	Water & Wastewater	Plastic or Glass	500 mL	Cool 6 C	48 Hours
Temperature	Water & Wastewater	Plastic or Glass	100 mL	None	Analyze Immediately
Turbidity	Water & Wastewater	Plastic or Glass	250 mL	Cool 6 C	48 Hours
<b>Metals</b>					
Metals, General	Water & Wastewater	Plastic or Glass, Rinsed w/ 1:1 HNO <sub>3</sub>	500 mL	HNO <sub>3</sub> to pH < 2	6 Months
	Solid	Plastic	50 mL	Cool 6 C	6 Months
Mercury	Water & Wastewater	Plastic or Glass, Rinsed w/ 1:1 HNO <sub>3</sub>	500 mL	1:1 HNO <sub>3</sub> to pH < 2	28 Days
	Solid	Plastic	50 mL	Cool 6 C	28 Days
Lead	Water & Wastewater	Plastic or Glass	1000 mL	HNO <sub>3</sub> to pH <2	6 Months
	Drinking Water	Plastic or Glass	1000 mL	HNO <sub>3</sub> to pH <2	6 Months

# COPY

Figure 3: continued

Parameter	Matrix	Container	Minimum Sample Size	Preservative	Holding Time
<b>Organic Compounds</b>					
GRO (Gasoline Range Organics), PVOC	Water	(4) 40 mL amber, glass bottle w/ 0.5mL 1:1 HCL, no headspace	40 mL	4 drops conc. HCL, Cool 6 C	14 Days
	Soil	(2) 60mL amber, glass bottle w/ 25mL MeOH	20- 35 grams	25 mL Methanol, Cool 6 C	21 Days
DRO (Diesel Range Organics)	Water	(2) 1 L Amber, glass bottle w/5.0 mL 1:1 HCL	1 L	5 mL 1:1 HCL, Cool 6 C	7 Days
	Soil	(2) 60 mL amber, glass bottle w/tare weight, no preservative	20- 35 grams	Cool 6 C	10 Days
BTEX	Water	(3) 40 mL glass vial w/ septa cap, zero headspace	40 mL	4 drops conc. HCL, Cool 6 C	14 Days
	Soil	(2) 4 oz Amber Glass w/ septa cap, zero headspace	25 grams	Cool 6 C	14 Days
VOCs:	Water & Wastewater	(4) 40 mL glass vial w/ septa cap, zero headspace	40 mL	0.5 mL 1:1 HCL	14 Days
	Soil	(2) 60 mL amber	20-35 g	25 mL MeOH	14 days
Pesticides, General	Water	(2) Amber Glass w/ TFE septa cap	1 Liter	Cool 6 C	7 Days
	Soil	(1) 4 oz Amber Glass	50 mL	Cool 6 C	14 Days
TCLP, Toxicity Characteristic Leaching	Water	1 L Amber Glass	1 Liter	Cool 6 C	14 Days
	Soil	(2) 4 oz Amber Glass	200 grams	Cool 6 C	14days

Solids & Semi-Solid - Microbiology	Matrix	Container	Minimum Sample Size	Preservative	Holding Time
Coliform, Fecal	Semi-Solid/Solid	Sterile Container	100 mL	Cool 6 C	8 Hrs EPA compliance 24 Hours

Solids & Semi-Solid - Inorganics and Nutrients					
% Solids	Semi-Solid/Solid	Plastic or Glass	100 grams	Cool 6 C	7 days

# COPY